



Serological Pipettes

1 ml – 70 ml



Product range of TPP serological pipettes.










An ideal combination is the use of TPP tissue culture flasks with TPP serological pipettes: you will experience total access to the corners of even the smallest flasks with our large pipette.



High capacities are obtained with reservoirs.

- Complete product range of serological pipettes from TPP: accurate calibrated volume from 1 ml – 70 ml.
- Single wrapped, sterile pipettes in stackable dispensers with dual opening mechanism:
 - individual dispense convenient through the tear open flab.
 - withdrawal of large quantities through the easy to open flip top.
- 94525: volume 25 ml + reservoir 15 ml = total capacity 40 ml
- 94550: volume 50 ml + reservoir 20 ml = total capacity 70 ml
- Nominal volume is significantly increased through reservoirs.
- Specially wrought tips of pipettes without burrs cast in one mould guarantee an excellent drip free handling.
- The smooth inside of the tubing without transition from top to tip prevents accumulation of residues. This provides further security in the routine.
- Slender form with a small outside diameter and a short working length! This optimises the ergonomically pipetting under the laminar flow hood.
- A low endotoxin level of < 0.1 EU/ml and a crystal clear transparency.
- An optimised mouthpiece with all pipettes conserves the pipette holder of your pipette controller.

For further technical information please visit us at www.tpp.ch.

Type	Product No.	Volume	Reservoir	Colour Code	Dimensions	Volume Graduations	Qty / Case
	94001	1 ml		yellow	Ø 4.8 x 270 mm	1/100 ml	400
	94002	2 ml		green	Ø 5.5 x 270 mm	1/100 ml	300
	94005	5 ml		blue	Ø 9.5 x 295 mm	1/10 ml	200
	94010	10 ml		orange	Ø 11.0 x 295 mm	1/10 ml	200
	94024	25 ml		red	Ø 16.0 x 300 mm	2/10 ml	100
	94525	25 ml	15 ml	red	Ø 13.5 x 345 mm	2/10 ml	60
	94550	50 ml	20 ml	purple	Ø 18.5 x 345 mm	1/10 ml	50

Your local dealer