





Tissue Culture Dishes

Large Surface and Secure Grip



Tissue culture dishes from TPP with crystal clear surfaces

Clock-numbering system inlaid in the base



Convenient peel-off sleeve

Proliferation

- A completely flat growth surface provides even and consistent cell growth resulting in optimal cell proliferation.
- Large growth surface, up to 10 % larger than available alternatives*.


Security

- Gripping ring: often imitated, never matched! The spikes of the protruding edge on the dish base ensure a secure grip of the dishes. Accidental lift-off of the lid becomes a thing of the past!
- The convenient peel-off sleeve enables the dishes to be withdrawn easily without compromising sterility.
- Secure and stable stacking of several dishes is facilitated by a raised stacking ridge on the lids and bases edges. Stacked dishes can be securely transported from A to B.
- As an ideal product combination for a convenient tissue harvest we recommend the use of the cell spatula or the cell scraper with twistable blade.

Documentation




- Numerical orientation by the clock-numbering system inlaid in the base: divides the dishes into quadrants for quick and easy orientation.
- Excellent microscopy through crystal clear surface, also ideal for microphotography.
- Yellow inscription field on the lid and frosted inscription field on the base provide convenient labelling.

For further technical information please visit us at www.tpp.ch.



Comparison of growth surfaces
 (*Source: internet)

Type	Product No.	Growth Surface	Internal-Ø	Dimensions	Material	Qty / Bag	Qty / Case
○	93040	9.2 cm²	34 mm	Ø 40 x 11 mm	PS	20	900
○	93060	22.1 cm²	53 mm	Ø 60 x 16 mm	PS	14	840
○	93100	60.1 cm²	87 mm	Ø 96 x 21 mm	PS	10	240
○	93150	147.8 cm²	137 mm	Ø 146 x 21 mm	PS	5	100

Type	Product No.	Length	Version	Width of Blade	Material	Qty / Bag	Qty / Case
	99002	240 mm	twistable	13 mm	PP/PE	1	150
	99003	300 mm	twistable	20 mm	PP/PE	1	150
	99010	195 mm	-	14 mm	PP	1	150

Your local dealer