



# BIOCHROM AG

# Neu von Biochrom

## A new exterior: All products now have new labels

[09-09-28] We have changed the appearance of our labels. The most important criterion: A colour guidance system helps you achieve the correct storage in your laboratory. Labels with red, green and blue marking each represent a temperature range.

**BIOCHROM AG**

**RPMI 1640 MEDIUM (1x)**  
 w 5.5 g/l NaCl    w 5 mg/l Phenol red    w 2.0 g/l NaHCO<sub>3</sub>  
 w 25 mM HEPES    w stable glutamine  
 LE ( Low Endotoxin )  
 Exp. 02 / 2009

**STORAGE**  
**+2-+8°C**

Cat. no.    **FG 1385**  
 Lot. No.    **0819 L**  
 Quantity    **500 ml**  
 For in vitro use only

ISO 9001 certified  
 Liquids: Solubilized in WFI  
 Made in Germany

Leonorenstr. 2-6 · D - 12247 Berlin



**STORAGE**  
**RT**  
 room temperature

**STORAGE**  
**+2-+8°C**

**STORAGE**  
**-20°C**

product to be stored at room temperature: normally at 15-25 °C	product to be stored in the cold room: in the specified temperature range	product to be stored deep-frozen: normally at ≤-20 °C +5 °C
---	--	--



### BIOCHROM AG

Leonorenstraße 2-6  
D- 12247 Berlin

Tel.: +49(0)30.779 90 60  
 Fax: +49(0)30.771 00 12  
 E-Mail: info@biochrom.de  
 Web: www.biochrom.de