



# BIOCHROM<sup>AG</sup>

# New from Biochrom


## Mycoplasma detection with low error rate and reduced effort: The new PCR-Kit Venor<sup>®</sup> GeM Advance

### Convenient Mycoplasma Detection

Venor<sup>®</sup> GeM Advance is a ready-to-go version of Venor<sup>®</sup> GeM. The reaction tubes included with the kit are pre-coated with appropriate dNTPs and primers to reduce the total assay time compared to the general protocol that requires individual loading of reaction tubes.

For additional convenience the gel loading buffer and dye are already included in the reaction tubes. After thermal cycling, the PCR can be loaded directly on the agarose gel.

**50 % introductory rebate  
only until 2009-31-12**



Cat. No.	W 11-7024	W 11-7048	W 11-7096	W 11-7240
Tests	24	48	96	240
List price	169,00 €	299,00 €	538,00 €	1139,00 €

### All benefits of Venor<sup>®</sup> GeM Advance at a glance:

- **Sensitive**  
Detection of 5 copies/μl of mycoplasma DNA that corresponds to 10 mycoplasma genomes per sample volume
- **Specific**  
Detection of more than 79 mycoplasma species
- **Reliable**  
with internal amplification control, optimized reproducibility by pre-aliquoted reagents.
- **Simple handling**  
PCR tubes included; just add rehydration buffer, polymerase, and sample.
- **Clear**  
A clear and easily interpretable result is obtained with one PCR reaction. Positive control included.

### Avoid mycoplasma contamination from the outset!

Zell Shield<sup>™</sup> protects cell cultures from contamination with bacteria, fungi and mycoplasma. Zell Shield<sup>™</sup> is a combination of unique and extremely powerful antibiotics. Common antibiotics are no longer required!

Cat. No.	W 13-0050	W 13-0150
Unit	50 ml	3 x 50 ml
List price	69,90 €	189,00 €



### BIOCHROM AG

Leonorenstraße 2-6  
D- 12247 Berlin

Tel.: +49(0)30.779 90 60  
Fax: +49(0)30.771 00 12  
E-Mail: [info@biochrom.de](mailto:info@biochrom.de)  
Web: [www.biochrom.de](http://www.biochrom.de)